

INSTRUCTIONS:

GERI-CHAIR - ARMCHAIR

1. Place product (A) in back of chair with bottom touching seat cushion and back loop strips (D), facing out.
2. Buckle behind chair and tighten (see Illustration 4). For a straight back chair, straps may be crossed if necessary (see illustration 2).
3. Place fleece (E) on front of cushion (A), and attach to Hook tabs (F) on back of cushion.
4. Attach (B) to (D) and press firmly.
5. Place (C) in the lap, one on each side of user, for arm and wrist support.

*Note (B) and (C) may be interchanged to meet the user's needs.

*To prevent skin irritation and discomfort, ALWAYS keep loop strips covered by fleece slip cover, on cushions not attached.




Illustration 2



Illustration 2 (straps may be crossed or pinned)

ALCO
Sales & Service Co.

**WHEELCHAIR**

1. Turn cushion (A) around so that the top is resting on the seat of the wheelchair, with the long back loop strips (D), on front facing out.
 2. Place the fleece (E) on front of cushion.
 3. Fold the top of cushion (A) under and place against chair back (see illustration 3).
 4. Bring straps behind chair, buckle and tighten (see illustration 4).
 5. For head and/or neck support, attach (B) across top to (A), adjust to comfort, press firmly.
 6. For shoulder support, attach (B) to (D) and press firmly.
- Use (B) or (C) in lap, either across knees, or under the arm and wrist, for support.

*To prevent skin irritation and discomfort, ALWAYS keep loop strips covered by fleece slip cover, on cushions not attached.



Illustration 3



Illustration 4