ALCO®
Sales & Service Co.

The Clear Solution To
Stop The Splash

INSTRUCTIONS FOR USE

CLEANING & MAINTENANCE:

- Looey® may be cleaned and disinfected using facility-approved disinfectant wipes or bleach.
- The frequency of cleaning and disinfection should follow your facility's established cleaning and disinfection protocols.

USE:

- Lower the Looey® shield and bedpan washing diverter arm to the down position.
- Pro Tip: Hold the vessel close to the shield and spray nozzle to minimize outward spray.
- Keep the faucet and Looey® in the upright position when not in use.
- Follow your facility's protocols for cleaning and disinfecting Looey® and clinical measuring containers.

WARNINGS:

- Do not sit, lean, pull, twist, or place objects on Looey®.
- The Looey® shield is not weight-bearing and should not be used to support any weight.
- Do not use Looey® if the shield is cracked or damaged, hardware is missing, or the shield is not securely attached to the diverter arm.
- Improper installation or use may result in inadequate splash protection.

www.alcosales.com | 800-323-4282 | 630-655-1900 | info@alcosales.com