

Purpose


Econo-Gel Pads are designed to relieve pressure and increase comfort for wheelchair users. High-density foam layers on top and bottom provide stability and cushioning. The gel-filled center inserts add pressure-relieving flotation.

Important Notice

Econo-Gel Pads have safety ties that must be secured to the wheelchair to prevent the cushion from slipping off the wheelchair seat.


Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser or mild disinfectant. The polyester cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.



Econo-Gel Pads - Standard Series

AL #	Cover	Size
AL-62360	Vinyl w/Polyester	16" x 16" x 2"
AL-62361	Vinyl w/Polyester	18" x 16" x 2"
AL-62362	Vinyl Only	16" x 16" x 2"
AL-62363	Vinyl Only	18" x 16" x 2"



Phone: 800-323-4282
Fax: 800-950-1167



Like us on Facebook at:
[facebook.com/ALCOSalesServiceCo](https://www.facebook.com/ALCOSalesServiceCo)



Follow us on Twitter at:
twitter.com/ALCOSales

info@alcosales.com
www.alcosales.com