

The advanced positioning of the AL-78405 allows the bed to go into manerous positions for

- Positively influencing clinically driven outcomes through multiple positioning options.
- Accommodating post-surgical reliab positioning to facilitate resident comfort
- Reduces caregiver strain and potential injuries to caregiver and residents caused by manually repositioning residents



Safety inspired embedded staff control in isothand • Angled for better visibility.

Angeot for before vanbeary.
Tamatasis locksets - backnet of semilant.

AL-78405

- use and lockout of embedded stall control. • Oaick lockout battons lock multiple functions
- simultaneously.
- leature (unless disabled), Patents Pending.



- Fail and Reverse Trendetenbarg positioning from 0 to 14 degrees, ideal for relation and subarute, provides therapeutic benefits that can aid in circulation and can body discomfort. Trendetenbarg also promotes gravity assist resident self positioning.
- Multiple positioning such as lossage chair and recliner style loss tissuity provides conduct and relaxation to the resident and also along for positioning of resident closer to the caregiver for proper care.
- Bed automatically passes when it returns to hevel from Trendeferdrang, Reverse Trendeferdran or chair position.
- Auto CPR functionality allows for ideal CPR positioning in emergency situations.



The Conflue Chair's privary's designed for resolutions with hourt and respiratory problems to whome a sumf-serce patients brings great relief. This chair can be made to serve as a fully reclised that we needed, hus can also be adjusted to form a chair. This "cardiac position" height to elevate the patient from the work equ which provides noted to the hings, personas circulations, and assists the gatent in coming from a fully process position to a string position without caudy harm or under strata.

Silversloutions¹⁴ within the powder-conting provide a untilaterobial surface protection 1018 99 5% of odes and stain causing bacteria Luchs for the expected Bi of the bed Nonalegenic and safe for the ervisionsent 104ph concertants on diversions given high efficacy 10ephen no special handling 1.cont effective solution to both in to last

800-323-4282 www.alcosales.com



Sealed ball bearings ensure maintenance-free, continuous, silent Designed for comfort and efficient operation

Height travel range: 7" to 30" Straight vertical rise protects walls and maintains footprint.

> Strong, no-sag grid deck sleep surface with sturdy cross braces. Safe work load: 600 pounds*



Synchronized high-low motors keep the bed level for proper therapeut positioning. When the head position is lower than the heart, lung volume may be compromised, causing airway problems, making it more difficult to breathe. Unlevel sleep surfaces may cause a risk of brain swelling and reduced systolic and diastolic pressures.



operation.

Intuitive two-pedal locking system is centrally positioned and easily visible from a distance to ensure that the bed is locked. Pedal locks are readily accessible at ALL bed heights and widths, even in the lowest caregiver back strain.





Grounded hospital-grade plug is a standard feature.



Bed provides 360° mobility at ANY height.

Directional mechanism on front caster is accessible at all heights for ease in rolling the bed. Timesaver Technology™ - fast motors save time.





Bed easily adjusts for use at 76° or 80° length. • Adjustable wall saver positions (76°, 80°, and for use with a trapeze) provide the smallest footprint at each length, while protecting the wall (included with bed).

- Adjustable mattress retainer securely supports mattress position for 76" or 80" length (included with bed).
- Edema foot ratchet to elevate the foot end of the bed (included with bed).



Maximum underbed clearance for use with overbed tables and lifts.

800-323-4282 www.alcosales.com

Ontional Half Length Assist Device and

Assist Bars provide a sturdy and secure hand hold to assist residents in and out of bed. Assist devices are removable without tools and meet FDA Entrapment Guidelines.



Warm-to-the-touch Half Length Assist Device AL-79178





Pivot Assist Bar

(tool-less) (bolt-on)

Pendant Hand Control and Holster

Easy operation contoured pendant is shock-resistant, splash-proof and durable. Pendant connects to either side of the bed to adapt to varying room layouts. Holster conveniently secures to bedrails and assist devices for easy access. Backlit hand control with underbed light is optional.

Battery Backup

Provides the ability to perform bed functions in the event of power outage. AL-69377 Battery backup AL-69378 Charger



sleep deck to 39" or 42" width for greater resident comfort. (Head and foot boards and staff control electronics not included.)



May be attached with bolts or quick release pins.



Warm-to-the-touch half-length assist device and

Move out of the way for lateral transfers.

 Pivot away from egress area. Provide toe clearance even in lowest

position to prevent staff injury. Hand pendant mounts facing inward or outward. and away from hand grab area.

nivot assist hars:

(tool-less) (bolt-on)

Extension Kit Transform your bed to 84" with

this optional 4-inch sleep deck extension. Standard 35" width for 39" or 42" widths



Trapeze Adapter

Optional trapeze adapter allows for easy attachment of trapeze system to headboard within easy reach for residents for repositioning themselves.



Trapeze Bar)

Wide Deck Kit

Optional wide kit adjusts the bed



www.alcosales.com

AL-78405

10000

Series Specifications

U.S. Patents - Multiple Patents Pending

Height Travel Range	7" to 30"
Safe Work Load*	
Mattress Deck Length	76"/80"
Mattress Deck Width	35"
Overall Bed Length (76"/80")	81.75"/85.75" with wall saver
Overall Bed Width	
Overall Bed Width	40.75" with half length assist devices
Maximum back angle	70°
Maximum knee angle	30°

Tool-less assembly minimizes set-up time.

Proudly Made in America



NEW Improved Warranty

- · 15 years on frame
- · 5 years on motors and control box
- · 3 years on other electronics

The Zenith9000 bed and accessories are in full compliance with the FDA guidelines for reducing the risk of bed entrapment, "Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment". Details can be found at: www.fda.gov.

For more information on the Zenith Series, please visit us online at www.grahamfield.com or call Basic American today.

* including resident, support surface, head/foot boards, and accessories.

Highest Safety and Quality Standards <u>Unrestricted</u> UL Listing

Built to conform to UL 60601-1, IEC601-2-38, and CAN/CSA C22.2 No. 601.1 the Zenith Series offers the following, ensuring the utmost resident and staff safety:

- · No pinch points within 200mm of bed perimeter
- Bed can be used as a transport device because it passes threshold clearance at ALL heights by 20mm
- · Splash-proof electronics
- · Load capacity and tip-stability tested and approved



800-323-4282 www.alcosales.com