# Powerheart<sup>®</sup> G5 + ICPR



The easy-to-use Intellisense<sup>®</sup> CPR Feedback (ICPR) device works in combination with your Powerheart G5 AED to delivery potentially life-saving CPR and defibrillation therapy.

## AL-88671

## About the Powerheart G5 + ICPR

Designed to be used by anyone, the Powerheart G5 AED with ICPR Feedback provides the guidance and shock power needed when someone goes into sudden cardiac arrest.

With prompts you can both see and hear, Intellisense CPR Feedback is the only CPR aid of its kind to combine compression rate and depth corrections to the 2015 AHA and ERC CPR Guidelines.



### **Powerheart G5 Features**



#### Intellisense CPR Feedback

Only ICPR delivers four distinct corrective prompts for rate, depth, and recoil



#### Always Rescue Ready®

Daily self-tests ensure that your device is ready to use



#### Rugged & Easy To Use

From novice to professional, it tells you exactly what to do



#### Dual Language Capibility

Switch from English to Spanish (voice and onscreen text prompts) with the push of a button!



#### Variable Escalating Energy

Subsequent shocks increase to a higher dose of energy

## Effective CPR Now in Your Hands AL-88671

## Introducing Intellisense<sup>™</sup> CPR Feedback

Performing effective CPR is crucial during a cardiac emergency and often leads rescuers to wonder, "Am I doing CPR properly?" With Intellisense CPR Feedback, providing effective CPR compressions is now within reach.

Thanks to a compact design, Intellisense CPR Feedback may be used on patients of all chest sizes. Motion sensor technology — the same technology used on most smart phones — provides rescuers with immediate voice prompts if their CPR compression rate and depth are outside the 2015 AHA-acceptable parameters.\* The easy-to-use, intelligent CPR device works in combination with your Powerheart® G5 automated external defibrillator to deliver potentially life-saving CPR and defibrillator therapy.

Discover today why effective CPR is now made easy.



#### Intelligent Technology, Effective CPR

- Each device features Rescue Ready<sup>®</sup> technology, embedded smart circuitry, and a motion sensor to self-test for presence and functionality.
- With one touch, CPR Feedback voice and text prompts are provided in an alternate language.
- Rescuer CPR performance can be captured on the rescue report.
- Rescue Ready self-testing data captured in AED Manager.



Supports AHA / ERC 2015

#### Simple to Prep, Simple to Use

- Experienced and first-time AED responders can perform effective CPR based on 2015 AHA/ERC guidelines.
- Simple place-and-press functionality requires no complicated preparation, adjustment, or delay.
- Non-slip surface and ergonomic shape ensures secure placement during CPR.
- Voice and text feedback provides corrective rate and depth instructions to help ensure effective CPR.
- If a rescuer does not want to use CPR Feedback, he or she may set it aside without it interfering with the rescue.



#### "Press Harder and Fully Release"

- Industry-first recoil analysis that provides voice and text feedback to indicate when compression depth needs adjustment. "Press harder and fully release" if compressions are shallow; and an emphasis on recoil. "Press softer" if compressions are too hard.
- With Intellisense CPR Feedback, the Powerheart G5 is the first AED to provide dual-language corrective prompting for both compression rate and depth in accordance with the 2015 AHA / ERC CPR Guidelines.

\*Based on European Resuscitation Council (ERC) and American Heart Association (AHA) 2015 Guidelines for Cardiopulmonary Resuscitation (CPR) and Emergency Cardiovascular Care (ECC) released on Oct. 15, 2015. \*\*For use on cardiac arrest patients eight years of age or older, or who weigh more than 55 lbs (25 kg).

ALCO Sales & Service Co. • 6851 High Grove Blvd., Burr Ridge, IL 60527 USA

800.323.47282 • Fax: 800.950.1167 www.alcosales.com • info@alcosales.com

## Visit www.alcosales.com

