


## Gel-Foam Heel Cushion

<u>AL #</u>	<u>Description</u>	<u>Unit</u>
□AL-62431	Gel-Foam Heel Cushion	Pair



### **An improved design for pressure sore prevention.**

The Gel-Foam Heel Cushion combines comfortable, resilient foam with gel to relieve heel pressure. The water-based gel insert cools to reduce skin damaging perspiration. Covered with a brushed polyester fabric that prevents skin shearing and friction, the Gel-Foam Heel Cushion is durable and completely launderable.

### **Easy to Use**

Place on patient's foot making certain that gel pad is positioned beneath heel. Secure Velcro<sup>®</sup> closure.

### **Easy Care**

Remove gel insert before laundering. (See illustration.) Wash and dry at temperatures that do not exceed 180°F (82°C). Do not use bleach.

