



Heel-Float
sold separately
AL-82941

AL#	Description	Size	Units
AL-82940	Bilateral Anti-Rotation Gripper™	11 3/8" x 8" x 3 1/2"	Each
AL-62432	Plastic Bilateral Anti-Rotation	Small/Medium	Each

Purpose

In the event it is advantageous to maintain the leg in a neutral position when lying in bed this bilateral anti-rotation gripper is recommended to be used in conjunction with most heel boots.

Easy to Use

First place the appropriate heel boot on the patient, then place the foot and boot into the Bilateral Anti-Rotation Gripper while lying on the bed. The Bilateral Anti-Rotation Gripper will release from the boot when the boot is lifted out of the Bilateral Anti-Rotation Gripper.

Easy Care

The Bilateral Anti-Rotation Gripper can be wiped clean with soap and water. If disinfected, avoid using chlorine (bleach) chemicals. Do Not Launder.

This product is guaranteed for 6 months.