

Re-Order Number	Description	Size
Call	Gel-Foam - Soft Base	16" x 16"
AL-84690-S	Foam Only - Soft Base	18" x 16"
Call	Foam only, Soft Base	20" x 16"
AL-84704-S	With Gel Pad, Soft Base	18" x 16"
Call	Gel-Foam - Soft Base	20" x 16"
AL-84690-F	Foam only, Firm Base	18" x 16"
Call	Foam only, Firm Base	20" x 16"
AL-84704-F	With Gel Pad, Firm Base	18" x 16"
Call	Gel-Foam, Firm Base	20" x 16"
AL-83937	Gel-Foam Cushion for Geri-Chair	20" x 18"



Purpose

Provides restraint-free positioning for wheelchair users with forward thrust and sacral sitting problems. Maintains pelvic neutrality. Designed as a comfortable alternative to pommel-style cushions, especially for males.

Important Notice

This cushion has tie straps and should be used at all times. Although this cushion is designed to help prevent the resident from sliding off a chair, the resident should be monitored as needed. If the resident should be monitored as

Easy Care

The outer cover can be wiped clean and is removable. It can be hand or machine washed, air or machine dried at temperatures less than 140°F (60°C). The inner cover is heat sealed and can be wiped clean with a standard disinfectants.

