




Contoured base reduces sling-seat hammocking. It rests below frame so user sits lower in chair

Position Plus[™] Contour Cushions

AL #	Size	Thrust Barrier Height	Unit
AL-84699-16	16"W x 16"D x 2.5"H	1	Each
AL-84699-18	18"W x 16"D x 2.5"H		Each
AL-84699-20	20"W x 16"D x 2.5"H		Each



Contoured wedge base controls forward sliding and reduces sling-seat hammocking. It rests below seat frame so user sits lower in chair

Position Plus[™] Contour Wedges


AL #	Size	Thrust Barrier Height	Unit
AL-84699-16W	16"W x 16"D x 3.5"H x 1"H	1	Each
AL-84699-18W	18"W x 16"D x 3.5"H x 1"H		Each
AL-84699-20W	20"W x 16"D x 3.5"H x 1"H		Each

Safety Strap keeps the cushion in place

Water-based gel relieves pressure and cools the skin to prevent breakdown caused by perspiration.

Side gel pods improves lateral stability by controlling left-right leaning.

High-density foam adds comfort and prevents bottoming out.



Rear gel pod provides pressure relief in the coccyx area.

Low-Shear I outer cover reduces skin-damaging friction.

Incontinent-proof vinyl cover prevents urine from entering cushion. Wipes clean.

Front gel pod controls forward sliding by eliminating posterior pelvic rotation.

Purpose

Foreward sliding and lateral leaning must be controlled to achieve proper positioning. Cushions and Wedges are designed with four gel-filled pods that act independently to keep the wheelchair user centered on the seat. These comfortable cushions and wedges reduce interface pressure, eliminate sling-seat hammocking, and prevent skin-damaging heat build-up.

Easy Care

The vinyl inner cover may be cleaned with a damp cloth. Spray-type cleaners and mild disinfectants may be used. Do not launder or dry cloth cover at temperatures below 180°F (82°C). If there isn't a gap, it may not be possible to apply the belt.

Caution

The Cushions and Wedges have safety straps that must be secured to the wheelchair to prevent them from sliding off the seat.