

A close-up, high-resolution photograph of a woman's face, showing her eyes, nose, and mouth. She has dark hair and is wearing a small diamond earring. The background is a soft, out-of-focus grey.

ALCO[®]
Sales & Service Co.



Healthcare Beds

Introduction

We produce and provide lifting solutions for people with special needs. Our products are made to serve as a natural part of life to all our users.

Human Care has a proud history in patient care. Since being founded in 1993, our twenty years of experience is expressed through a broad selection of premium products offering high customer value and satisfying all customer needs.

Our main product areas are: Healthcare Beds, Mobility Aids, Lifting Solutions, Convertible Chairs, and Bathroom Safety Products for people with reduced mobility.

With headquarters in Stockholm, Sweden, and subsidiaries in the USA, Canada, Holland and Australia, we collaborate with distributors to provide solutions in more than twenty countries.

Our promise is safety, quality and usability. In all of our products, we exercise only the highest quality to make sure that every detail, every aspect, every function will deliver that promise made to our users.

Welcome to Human Care.

Distributed by

ALCO[®]
Sales & Service Co.



Content

Healthcare Beds	Page
Reduce fall injuries	4
Care and quality	5
FloorLine- <i>i</i> Plus	7
FloorLine- <i>i</i>	8
FloorLine LTC	9
Accessories	10

Your Life. Your Way.

That is what we have been offering people with reduced mobility for over 20 years. Your Life. Your Way. Our focus is simple. We work for our users. Our products serve to be a natural part of life to anyone who might need them. Meeting that need is what drives and motivates us. Our mission is to improve the quality of care for people with special needs. Being taken care of by someone else – to literally be in someone else's hands – is a special situation. It demands trust, respect, a fine touch and – not least – safety for both user and caregiver.

Our solutions provide reliability, high quality and innovation. With low-maintenance requirements and efficient installation procedures, we provide a state-of-the-art offering for both users and caregivers.

In every detail, in every aspect, our products deliver what they promise.

Floor-level-beds reduce fall injuries without the use of restraints*

10 cm (Four inches). That is the difference between “falling” – or “rolling” – out of bed! So small a distance, yet so large a difference between extended hospital stays and cost – or their avoidance. Our range of Healthcare beds aims to reduce fall injuries and provide the opportunity for excellent care.*

Dignity. Functioning at a floor level, restraints are minimised, thereby greatly enhancing the dignity of care. One may elect to call it ‘dignity’. We choose to simply call it ‘care’. Human Care.

Fall injuries cause pain – and cost

Patient falls are widely recognized as a significant healthcare risk, as well as a cost. Pain, cost, complication, extended rehabilitation stay – and even death – are just some of the dire consequences of a bed fall.

* Barker A., Kamar J., Tyndall T. and Hill K. (2013) Reducing serious fall-related injuries in acute hospitals: are low-low beds a critical success factor? Journal of Advanced Nursing 69 112-121

The risk of falling. Patients who are at risk of falling are not protected from the impact of a fall if the bed only descends to a height of 20 cm [8 inches]. That distance to the floor – although appearing small and harmless – can still lead to serious injuries.

Avoid fall injuries. On a floor-level-bed with a minimum height of less than 10 cm [4 inches], the gap from the bed to the floor is significantly lower, which reduces the force of the fall. Thus, the risk of bed-fall injuries is reduced, offering safety – as well as dignity – to the user. Our range of premium beds – featuring all state-of-the-art functions expected in today's leading healthcare equipment – reduces patient's risk of bed fall injuries to an unprecedented low.

Improve the working position. A challenge in today's healthcare is the working position of caregivers. To resolve this problem, the floor-level-beds may be raised to any ergonomically suitable height up to 799 mm [31 inches] from the floor, creating an optimum working height that reduces back injuries risks for all caregivers.



” Our mission is to improve the quality of care and freedom of movement for people with special needs ”

Care and quality

A user of a product from Human Care is a person who feels very confident with the product's quality and reliability. At all times. Every step of the way. The reason behind this is our profound commitment to safety, quality and usability.

Firstly, our products are designed and manufactured to meet the requirement of relevant ISO and EN standards. Secondly, every single product is given extensive attention at every stage of the design process. Thirdly, in developing our products, our design team collaborates with globally leading physiotherapists, creating state-of-the-art products which all share the same ultimate goal: the benefit and enjoyment of their users. Allowing them to live life in the best way possible. Their way.

” 10 cm (4 inches) makes the difference between falling or rolling out of bed – These 10 cm reduce the risk of extended hospital stays and direct costs. ”



FloorLine-i Plus, FloorLine-i and FloorLine LTC beds are all electric true floor-level-beds that go down to an ultra-low level less than 10 cm (4 inches) above the floor as measured from the top of the mattress platform.

A close-up photograph of an elderly man with wrinkled skin and white hair, lying in a hospital bed. He is wearing a light blue hospital gown. A clear IV tube with a blue connector is visible near his neck. The background shows white hospital linens.

Floor-level-beds

electric floor-level-beds help alleviate safety concerns and are also a convenient way of taking care of patients and residents in hospitals and nursing/long-term care facilities. Placed at a true floor level they provide a vital care factor for those who may be at risk of falling out of bed and when raised they provide an ergonomically safe working height for any caregiver.

AL-85306

FloorLine-i Plus

Floor-level-bed for hospital environment

- A complete management solution
- Maximum height adjustable
- Future proofed for the addition of accessories.

The aesthetically appealing FloorLine-i Plus is a complete care management solution. With its extensive range of features such as trendelenburg, cardiac chair position and electronic CPR function, this bed offers our highest level of care. With its open bus digital control box it is also future proofed for adding even more accessories.

The FloorLine-i Plus is the perfect solution for preventing fall injuries. Its minimum height from the floor is just 10 cm (3.9"). At the same time no compromises are made for excellent ergonomic working height reaching 80 cm (31.4").



Part No. (including hand control and end boards)	Platform Position min- max	Sleeping Dimensions L x W	Bed Dimensions L x W	Functions	Compliance Standards	Maximum User Weight/ Safe Working Load	Product Weight
AL-85306	10 - 80 cm (3.9" - 31.4")	2000mm x 900 mm (78.7" x 35.4") With Extension Kit: 2175mm x 900mm (85.5" x 35.5")	2325mm x 925mm (91.5" x 36.4") With Extension Kit 2500mm x 925mm (98.5" x 36.5")	Twin castors Battery Back Up, Auto Contour Button for Cardiac Chair Posi- tion, Handset w Lockout	IEC 60601-2- 52:2009 EN1970:2002 AS/NZS 3200.2.38:2007	Maximum user weight: 185 kg (408 lbs) Safe working load: 250 kg (551 lbs)	140 kg (309 lbs)

FloorLine-i Plus can be fitted with a variety of optional accessories such as dedicated side-rails, self-help pole and length extension kits.

AL-85305

FloorLine-i

Floor-level-bed for use in hospitals and long-term care facilities

- Minimum height from the floor of only 10 cm
- Can be used both in hospitals and long term care facilities
- Minimizes patient restraints and use of side rails

The FloorLine-i is a sound purchase. It offers a variety of functionality such as trendelenburg and cardiac chair position providing excellent care.

The FloorLine-i is an optimal floor-level-bed solution for use both in hospitals and long-term care facilities. Its minimum height from the floor is just 10 cm (3.9") and at the same time no compromises are made for caregivers' ergonomics.



Part No. (including hand control and end boards)	Platform Position min- max	Sleeping Dimensions L x W	Bed Dimensions L x W	Functions	Compliance Standards	Maximum User Weight/ Safe Working Load	Product Weight
AL-85305	10cm - 80cm (4" - 31.5")	2000mm x 900mm (78.7" x 35.5")	2325mm x 925mm (91.5" x 36.5")	Twin castors Battery Back Up, Auto Contour Button for Cardiac Chair Position, Hand-set w Lockout	IEC 60601-2-52:2009 EN1970:2002 AS/NZS 3200.2.38:2007	Maximum user weight: 185kg (407lbs) Safe working load: 250kg (550lbs)	138 kg (303 lbs)
		With Extension Kit 2175mm x 900mm (85.5" x 35.5")	With Extension Kit 2500mm x 925mm (98.5" x 36.5")				

FloorLine-i can be fitted with a variety of optional accessories such as dedicated side-rails, self-help pole and length extension kits.

AL-85307

FloorLine LTC

Floor-level-bed for long term care facilities

- Lowest height from the floor of only 9.5 cm (3.7")
- The optimal choice for long term care
- Easy to fold and store efficiently

The FloorLine LTC with its home environment design is the optimal choice for long term care facilities. Its one-piece design makes it easy to fold and store away efficiently.















The FloorLine LTC is the perfect solution for preventing fall injuries. Its minimum height from the floor is just 9.5 cm(3.7") providing a safe and considerate way of taking care of people who may be at risk of falling out of their bed.



Part No.	Platform Position min- max	Sleeping Dimensions L x W	Bed Dimensions L x W	Functions	Compliance Standards	Maximum User Weight/ Safe Working Load	Product Weight
AL-85307	9.5 cm - 70 cm (3.7" - 27.5")	2000mm x 900mm (78.7" x 35.5") With Extension Kit 2175mm x 900mm (85.5" x 35.5")	2150mm x 905mm (84.5" x 35.5") With Extension Kit 2325mm x 905mm (91.5" x 35.5")	Auto Contour Button for Cardiac Chair Position, Handset w Lockout {Batteribackup as accessory, part no: 75805}	IEC 60601-2- 52:2009 EN1970: 2002 AS/NZS 3200.2.38: 2007	Maximum user weight: 150kg (330lbs) Safe working load: 200kg (440lbs)	106 kg (234 lbs)

FloorLine LTC can be fitted with a variety of optional accessories such as dedicated side-rails, self-help pole, Back Up Battery and length extension kits.

Accessories

<p>Hanger for air mattress pump Attachment to hold air mattress pump. To be mounted on the footend of the bed.</p> <p>AL-85308</p>		<p>Side rails 1/2 side rail</p> <p>AL-85315</p>	
<p>Care Assist Rail Designed to help patients safely in and out of bed.</p> <p>AL-85309</p>		<p>Side rails Split side rail (1/2 + 1/4 length side rail)</p> <p>AL-85316</p>	
<p>Push handle The push handle can be mounted on the foot end of the bed to improve the operation of the bed.</p> <p>AL-85310 (Included in the product for FloorLine-i Plus, accessory for FloorLine-i)</p>		<p>Side rails 3/4 side rail</p> <p>AL-85317</p>	
<p>IV pole (adjustable in height) Adjustable IV pole that can be mounted on both ends of the bed.</p> <p>AL-85311</p>		<p>Bumper The bumper protects the head of the bed and are mounted on the lower support bracket between the wheels.</p> <p>AL-85318-02 (FloorLine LTC) AL-85318-01 (FloorLine-i and FloorLine-i Plus)</p>	
<p>Extension kit Extension kit that extends the bed, with 175 mm to 2175 mm to accommodate more people.</p> <p>AL-85312-02 (FloorLine LTC) AL-85312-01 (FloorLine-i and FloorLine-i Plus)</p>		<p>Swing-away trapeze The free-standing swing-away trapeze allows patient to reposition themselves. No tools needed for assembling.</p> <p>AL-85319</p>	
<p>Oxygen Holder An attachment for placing an oxygen bottle in C size.</p> <p>AL-85313</p>		<p>Bed Raising Block (60 pcs/ box) The bed raising block is used for a simple way to raise or lower the seating or lying surface on beds, headboards, chairs or tables. May be stacked maximum 4 on top of each other.</p> <p>AL-85320</p>	
<p>Self help pole The self-help pole helps patients to get up in the bed and change positions.</p> <p>AL-85314-02 (FloorLine LTC) AL-85314-01 (FloorLine i and FloorLine-i Plus)</p>		<p>Ring for Bed Raising Block The bed raising block can be complemented with the ring 895-1 for rounded and narrow legs (less than 35 mm) to obtain a solid filled block.</p> <p>AL-85321</p>	



The FloorLine LTC with its home environment design is the optimal choice for long-term care facilities

Distributor:



6851 High Grove Blvd.
Burr Ridge, IL 60527

info@alcosales.com
www.alcosales.com
800.323.4282

