

### AL-82372

## Chair Scale for Weighing While Seated



- The simple solution for weighing wheelchair patients.
- Fold-up armrests and footrests make sitting down easy.
- Brakes on wheels for secure stance.
- Intuitive operating elements simplify use.
- All essential functions for maximum usage.

#### Relief for patient and medical personnel.

tanding on a scale is a problem for physically disabled or long-term care patients. Chair scales are of great assistance to patients and medical staff in nursing homes, hospitals, geriatric and orthopedic departments. The AL-82373 offers a sound, user-oriented solution with all the functions that a chair scale should offer.

# rterm care edical staff nts. The ctions that

#### TARE and auto-HOLD.

Of course the AL-82373 is also equipped with the important TARE function which allows every newly added weight to be ascertained separately and the auto-HOLD function which lets medical personnel take care of the patient first and read out the results later. All these advantages plus well-known workmanship and durability make the scale a genuine product. It fulfills all the expectations people have of the world leader in medical weighing and measuring.

#### Everything that counts.

With a capacity of 440 lbs and a graduation of 0.2 lbs, the AL-82373 fulfills all basic medical standards. Added to that are the typical advantages such as fold-up armrests and foot-rests, secure wheel brakes, an ergonomically shaped seat and smooth running rollers, which, together with the ergonomically shaped handle, make it very easy to push the chair scale.

#### Simple operation.

Thanks to the fold-up armrests and footrests, elderly and heavy patients can take a seat on the scale easily and comfortably. It is just as easy for personnel to use the scale because the batteries eliminate the need for connection to an electrical outlet and keep on working through up to 4000 weighings. No time is lost to a learning curve – a user can operate the scale intuitively. Display ele-ments are conveniently located in the handle and can be set to show kilograms, pounds or stones.









To make sitting down and standing up easier, the armrests and footrests fold up out of the way.



The practically positioned display is easy to

#### AL-82373

#### **Technical Data**

- Capacity: 440 lbs / 200 kg
- Graduation: 0.2 lbs / 100 g
- Dimensions (WxHxD):
   22.17 x 35.31 x 37.64" / 563 x 897 x 956 mm
- Weight: 49.42 lbs / 22.4 kg
- Power supply: Power adapter
- Functions: TARE, auto-HOLD,
   lbs/kg switch-over, automatic switch-off

