


AL-90374-02 & AL-90374-04 - Noti-Fi™ Wireless Smart Device Text Notification with Reporting Package



Description

The AL-90374-02 & AL-90374-04 Noti-Fi text alerting package for Android smart devices is used with the Pro-Alert 570 Nurse Call System to route active resident and staff calls over the facility's LAN to wireless smart devices carried by caregivers. In addition, it provides archiving and analysis reports of resident and staff call and response activities.

The smart devices load the Noti-Fi application over the LAN from the PC running the Noti-Fi software. Caregiver's smart devices are assigned in Noti-Fi to receive Pro-Alert 570 call alerts individually or for caregiver teams to receive alerts simultaneously. Caregivers can selectively ACCEPT an active call which in turn notifies other caregivers on the team who is responding to the call. Overtime reminders automatically alert caregiver(s)' smart device when a call has not been cleared within a pre-set time period. From their smart device, caregivers manage when they receive call alerts based on if they are on or off duty and if they are on break.

The AL-90374 is a turnkey package including all required software preconfigured on an all-in-one PC/monitor, which includes a keyboard, mouse, and license key for all required software. From the Reporting PC, staff generate reports customized based on all rooms or specific rooms/beds, all day or specific times, a single day or range of dates, and call priorities or specific call priorities. The included reports are: Summarized Call Statistics, Hourly Call Statistics, Summarized Activity, Detailed Activity, Activity Exceptions, and System Wide Activity. The reports are available in PDF, CSV, and XLS formats. The AL-90374 includes licensing for two smart devices; additional licensing, in increments of five devices, is available through the AL-90375 supporting up to an additional fifty smart phone devices.

When ordering, specify AL-90374-02 for system of up to two AL-90379T Main Control Units or AL-90374-04 for system of up to four AL-90379T Main Control Units.

Features

- Receive text alerts of active patient and staff calls on Android*-based phones or tablets
- Utilizes the facility's LAN to transmit call activity; no Internet connection is required
- At caregiver's smart device, a full display of active calls including room, bed, call priority, and call duration
- Display of at least eight calls at once with automatic prioritization based on call priority and duration
- Supports team nursing with ability to ACCEPT or REJECT calls
- Active call alert with tone and/or vibrate of smart device
- Calls display on locked smart phone screen
- Logs all resident to staff assignments and interactions on the Pro-Alert 570 Nurse Call system for future report generation
- Includes a set of standard reports with the ability to adjust each report's parameters
- No PC software installation required – all software comes preloaded on an all-in-one PC

PRODUCT SUBJECT TO CHANGE WITHOUT NOTICE

* Android version 6.0 and later supported

AL-90374E-02 & AL-90374E-04 - Noti-Fi™ Wireless Smart Device Text Notification with Reporting & Email Output Package



FEATURES

- Receive text alerts of active patient and staff calls on Android*-based phones or tablets
- Utilizes the facility's LAN to transmit call activity; no Internet connection is required
- At caregiver's smart device, a full display of active calls including room, bed, call priority, and call duration
- Display of at least eight calls at once with automatic prioritization based on call priority and duration
- Supports team nursing with ability to ACCEPT or REJECT calls
- Active call alert with tone and/or vibrate of smart device
- Calls display on locked smart phone screen
- Logs all resident to staff assignments and interactions on the Pro-Alert 570 Nurse Call system for future report generation
- Includes a set of standard reports with the ability to adjust each report's parameters
- Includes automatic report generation and emailing option
- No PC software installation required – all software comes preloaded on an all-in-one PC

Description

The AL-90374E-02 & AL-90374E-04 Reporting Package with Noti-Fi text alerting package for Android smart devices is used with the Pro-Alert 570 Nurse Call System to route active resident and staff calls over the facility's LAN to wireless smart devices carried by caregivers. In addition, it provides archiving and analysis reports of resident and staff call and response activities.

The smart devices load the Noti-Fi application over the LAN from the PC running the Noti-Fi software. Caregiver's smart devices are assigned in Noti-Fi to receive Pro-Alert 570 call alerts individually or for caregiver teams to receive alerts simultaneously. Caregivers can selectively ACCEPT an active call which in turn notifies other caregivers on the team who is responding to the call. Overtime reminders automatically alert caregiver(s)' smart device when a call has not been cleared within a pre-set time period. From their smart device, caregivers manage when they receive call alerts based on if they are on or off duty and if they are on break.

The AL-90374E is a turnkey package including all required software preconfigured on an all-in-one PC/monitor, which includes a keyboard, mouse, and license key for all required software. From the Reporting PC, staff generate reports customized based on all rooms or specific rooms/beds, all day or specific times, a single day or range of dates, and call priorities or specific call priorities. The included reports are: Summarized Call Statistics, Hourly Call Statistics, Summarized Activity, Detailed Activity, Activity Exceptions, and System Wide Activity. The reports are available in PDF, CSV, and XLS formats. The email output functionality allows reports to be scheduled for future generation and automatically emailed to a list of recipients. The AL-90374E includes licensing for two smart devices; additional licensing, in increments of five devices, is available through the AL-90375 supporting up to an additional fifty smart phone devices.

When ordering, specify AL-90374E-02 for system of up to two AL-90379T Main Control Units or Model AL-90374E-04 for system of up to four AL-90379T Main Control Units.

PRODUCT SUBJECT TO CHANGE WITHOUT NOTICE

* Android version 6.0 and later supported