

## AL-80522

### Safety Auto-Reset Monitor

### Our Most Popular Fall Monitor



### Our Most Popular Monitor!

AL-80522

### SAFE w/ "Always ON" Technology

The AL-80522 Monitor is designed for use with all Smart Caregiver corded sensor pads, floor mats and seat belts. Simply plug a sensor pad into the monitor and place the pad under the resident in a chair or bed. When the resident gets up and pressure is removed, the alarm will sound notifying caregivers that the resident is on the move. This alarm features a Safety Auto-Reset™ Function which will automatically reset and silence the alarm when pressure is re-applied to the sensor pad.

### Optional TamperProof™ Setting

When set to TamperProof™, the Monitor alert can ONLY be reset using the provided Caregiver Key™. Selecting this setting ensures that only the designated caregiver can silence the alarm.



### Key Features:

- Full-Range Volume Control
- Nurse Call Capability - enables you to use with existing hard wired nurse call system
- Optional Setting: TamperProof™ Caregiver Key allows control of who is able to reset the monitor
- Batteries not included (3-AA)
- Optional AC adapter

[info@alcosales.com](mailto:info@alcosales.com)

800.323.4282

[www.alcosales.com](http://www.alcosales.com)

