

TECHNICAL SPECIFICATIONS

SINGLE PATIENT STATION


AL-84673

ALCO
Sales & Service Co.


*Call Placement Indicator
Call Reset Button
Cordset Removal Triggers Call Placement
Standard Cordset Receptacle
UL Listed*



Front



Rear



Side

DESCRIPTION

Wall recessed or surface mounted single patient station for use with remote call origination cordset. Removal of the cordset originates a call placement.

RECEPTACLE

One cordset receptacle wired to accept a single 1/4" cordset (normal call, geriatric, explosion proof) with nurse call origination facility.

CONTROLS AND INDICATORS

One combination white call placement indicator and call cancel button.

INSTALLATION

Color-coded, prewired terminations provided. Mount using only the two 6-32 x 1" nylon screws provided.

HOUSING AND FINISH

Faceplate: Flame retardant ABS plastic rated 94V-O by UL, finished in off-white supplied with nylon mounting screws.

W41 Backbox: metal backbox supplied by others.

DIMENSIONS

Faceplate: 2-3/4" (6.99cm) wide x 4-1/2" (11.43cm) high.

W41 Backbox: 4-11/16" (11.91cm) square electrical box with one-gang adapter for a total depth of 2-7/8" (7.30cm), supplied by others.

LISTINGS

The AL-84673 is UL and CUL listed Hospital Signaling Equipment / Nurse Call Equipment under File E157381.

ARCHITECTS' AND ENGINEERS' SPECIFICATIONS

AL-84673 SINGLE PATIENT STATION

The wall recessed or surface mounted single patient station shall be for use with remote call origination cordset.

The patient station shall be provided with a white combination call placement indicator and call cancel button. A receptacle shall be provided, to accept single 1/4" type cordsets. A call shall be originated if the cordset is removed from the receptacle.

The unit shall be provided with color-coded prewired terminations for ease of installation. The faceplate shall be flame retardant ABS plastic rated 94V-O by UL. Nylon Screws shall be provided to mount the faceplate.

The unit shall be UL and CUL listed Hospital Signaling Equipment / Nurse Call Equipment under File E157381.

The faceplate shall be finished in off-white. The W41 backbox shall be supplied by others. The dimensions shall be: **Faceplate:** 2-3/4" (6.99cm) wide x 4-1/2" (11.43cm) high, **W41 backbox:** 4-11/16" (11.91cm) square electrical box with one-gang adapter for total depth of 2-7/8" (7.30cm), supplied by others.

